Gebrauchsinformation

ABO-trol[®] forte II

NUR ZUR IN-VITRO-DIAGNOSTIK

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ABO-trol[®] forte II

FOR IN-VITRO-DIAGNOSTIC USE ONLY

50 bedside-testcards ready for use for verification of ABO-compatibility prior to blood transfusion. Pre-dropped fluid antisera with 70 µl of each monoclonal anti-A (mouse) and anti-B (mouse) sufficient for one determination.

Clone designation anti-A: celline BIRMA-1
anti-B: celline LB-2

Manufacturer:
ANTITOXIN GmbH
Industriestrasse 88,
69245 Bammental

CE 0483

Testsera anti-A and anti-B used for ABO-trol[®] forte II are produced from monoclonal antibodies (mouse) that are distinguished by quick and strong reactions also with A-subgroups and immature A and B antigens.

Introduction

„Right before a transfusion the transfusing physician himself or under his direct supervision the ABO-identity test (bedside-test) has to be done with the recipient (e.g. with testcard). The test serves to confirm the results for ABO-bloodgrouping antigens of the recipient obtained from previous tests. The ABO-identity test additionally can be carried out with the erythrocyte concentrate that should be transfused”. (literature: „Richtlinien zur Blutgruppenbestimmung und Blut-transfusion“, 2000).

Product

ABO-trol[®] forte II consists of a plastic card with 4 trough-shaped cavities acting as reaction chambers. Each of these chambers is filled with 70 µl of antiserum, blue coloured anti-A on the left and yellow coloured anti-B on the right. The reaction chambers filled with antisera are sealed with a transparent film. The chambers are intended for confirmation testing of patient erythrocytes and donor erythrocytes relating to ABO compatibility.

Performance of testing

For identity testing erythrocyte concentrate, native blood and blood with anticoagulatives can be used.

1. Collect blood from the patient
2. Penetrate the transparent film of chambers „anti-A“ and add 1 single small drop of patient erythrocytes or donor erythrocytes free falling into these chambers.
3. As described under point 2 add 1 single small drop of patient erythrocytes or donor erythrocytes free falling into chambers “anti-B”.

4. Mix erythrocytes and antisera by gentle striking.
5. Check for agglutination (= positive reaction) after 10 seconds.

Interpretation

ABO compatibility is ascertained when bloodgroup of determined patient erythrocytes is identical with bloodgroup on label of the bag for transfusion.

The following schedule shows compatibility of ABO-different erythrocytes with poor plasma.

Anti-A	Anti-B	Compatible donor erythrocytes
pos.	neg.	A or O
neg.	pos.	B or O
pos.	pos.	AB, A, B or O
neg.	neg.	O

A different determination of bloodgroup of patient erythrocytes by this agglutination method compared to the signalized bloodgroup on the label of the bag for transfusion means that the bloodgroup of the bag has no identity with the bloodgroup of recipient.

Limitations of the procedure

False positive reactions may occur by: bacterial contamination of serum and/or of erythrocytes, Wharton'sche Sulze, antibodies diverted against colour ingredients or preservatives, medicinal therapy or pathological discrepancies, dextrans or other plasma expanders.

False negative results may occur by: drying up of testserum in reaction chambers or false concentration of erythrocytes.

False positive or false negative results lead to discrepancies between bedside test and previous test results. In this case ABO-compatibility is not confirmed and further tests have to be carried out with different test systems.

Storage and stability

ABO-trol[®] forte II is to be stored at + 2 to 8 °C. Under adequate storage the antisera keep their reactivity up to shelf life printed on the card.

Remarks

As preservative the testsera in ABO-trol[®] forte II contain sodium azide (<1 g/l). Harmful if swallowed, avoid contact with skin and mucous membranes

Warning:

These reagents were not tested for murine viruses. Do not pipette with mouth.

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